


Quantitative Depth Profile (QDP) Analysis of 40 GB Hard Disk on Glass Substrate

GDS Applications Lab • LECO Corporation, Saint Joseph, Michigan USA



GDS-850A analysis showing three replicate analyses of coating layers on a laptop hard disk.

Parameters

Grimm-Style RF 4 mm Lamp
5 W, 1000 V
10 s, 1000/s

For questions on this analysis e-mail us at:
info@leco.com

For a complete listing of snapshots and performance notes visit us on the web at www.leco.com in the Spectroscopy section of the Applications Library.

